


Product Datasheet

Insulin protein (Active) (orb85999)

Description	Insulin protein (Active)
Species/Host	Yeast
Conjugation	Unconjugated
Tested Applications	WB
Target	Insulin
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at -20°C or - 80°C storage. Avoid repeated freeze / thaw cycles
Note	For research use only
Purity	98% , Affinity chromatography
MW	14.3 kDa
Expiration Date	6 months from date of receipt.



Western blot analysis of Insulin protein