

Product Datasheet

hoxb8b antibody (orb860083)

Description	Rabbit polyclonal antibody to hoxb8b. All of our unique antibodies are supplied with the following: 2...
Species/Host	Rabbit
Reactivity	Fish
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Danio rerio (Zebrafish) (Brachydanio rerio) hoxb8b protein
Target	hoxb8b
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8JH55
Expiration Date	12 months from date of receipt.