

Product Datasheet

POMGNT2 antibody (orb860801)

Description	Rabbit polyclonal antibody to POMGNT2. All of our unique antibodies are supplied with the following: ...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Homo sapiens (Human) POMGNT2 protein
Target	POMGNT2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8NAT1
Expiration Date	12 months from date of receipt.