

Product Datasheet

rab11fip4a antibody (orb861000)

Description	Rabbit polyclonal antibody to rab11fip4a. All of our unique antibodies are supplied with the following...
Species/Host	Rabbit
Reactivity	Fish
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Danio rerio (Zebrafish) (Brachydanio rerio) rab11fip4a protein
Target	rab11fip4a
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q3LGD4
Expiration Date	12 months from date of receipt.