

Product Datasheet

SCR antibody (orb861157)

Description	Rabbit polyclonal antibody to SCR
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Metridium senile (Brown sea anemone) (Frilled sea anemone) Small cysteine-rich protein
Target	Small cysteine-rich
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P0DL60
Expiration Date	12 months from date of receipt.