

Product Datasheet

SCR8 antibody (orb861192)

Description	Rabbit polyclonal antibody to SCR8
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Orbicella faveolata (Mountainous star coral) (Montastraea faveolata) Small cysteine-rich protein 8 protein
Target	Small cysteine-rich protein 8
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	B2ZG38
Expiration Date	12 months from date of receipt.