

Product Datasheet

SCR5 antibody (orb861226)

Description	Rabbit polyclonal antibody to SCR5
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Orbicella faveolata (Mountainous star coral) (Montastraea faveolata) Small cysteine-rich protein 5 protein
Target	Small cysteine-rich protein 5
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	C1KIZ4
Expiration Date	12 months from date of receipt.