

Product Datasheet

tbc1d8b antibody (orb861574)

Description	Rabbit polyclonal antibody to tbc1d8b. All of our unique antibodies are supplied with the following: ...
Species/Host	Rabbit
Reactivity	Fish
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Danio rerio (Zebrafish) (Brachydanio rerio) tbc1d8b protein
Target	tbc1d8b
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	BOROW9
Expiration Date	12 months from date of receipt.