

Product Datasheet

TBPL2 antibody (orb861622)

Description	Rabbit polyclonal antibody to TBPL2.All of our unique antibodies are supplied with the following: 20...
Species/Host	Rabbit
Reactivity	Frog
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Xenopus laevis (African clawed frog) tbpl2 protein
Target	tbpl2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q5UE94
Expiration Date	12 months from date of receipt.