

Product Datasheet

tnfaip8l2a antibody (orb861655)

Description	Rabbit polyclonal antibody to tnfaip8l2a. All of our unique antibodies are supplied with the following...
Species/Host	Rabbit
Reactivity	Fish
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Danio rerio (Zebrafish) (Brachydanio rerio) tnfaip8l2a protein
Target	tnfaip8l2a
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q1ECV8
Expiration Date	12 months from date of receipt.