

Product Datasheet

UTS2B antibody (orb861835)

Description	Rabbit polyclonal antibody to UTS2B. All of our unique antibodies are supplied with the following: 20...
Species/Host	Rabbit
Reactivity	Fish
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Danio rerio (Zebrafish) (Brachydanio rerio) uts2b protein
Target	uts2b
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q7ZZY8
Expiration Date	12 months from date of receipt.