

Product Datasheet

atp6v0d1 antibody (orb861838)

Description	Rabbit polyclonal antibody to atp6v0d1. All of our unique antibodies are supplied with the following:...
Species/Host	Rabbit
Reactivity	Fish
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Danio rerio (Zebrafish) (Brachydanio rerio) atp6v0d1 protein
Target	atp6v0d1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q6PGV1
Expiration Date	12 months from date of receipt.