

Product Datasheet

TGIF2LX antibody (orb862170)

Description	Rabbit polyclonal antibody to TGIF2LX. All of our unique antibodies are supplied with the following: ...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Homo sapiens TGIF2LX protein
Target	TGIF2LX
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8IUE1
Expiration Date	12 months from date of receipt.