

Product Datasheet

At3g57050 antibody (orb862240)

Description	Rabbit polyclonal antibody to At3g57050. All of our unique antibodies are supplied with the following...
Species/Host	Rabbit
Reactivity	A. thaliana
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Arabidopsis thaliana At3g57050 protein
Target	At3g57050
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P53780
Expiration Date	12 months from date of receipt.