

Product Datasheet

Major allergen Can f 1 antibody (orb862255)

Description	Rabbit polyclonal antibody to Major allergen Can f 1. All of our unique antibodies are supplied with ...
Species/Host	Rabbit
Reactivity	Canine
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Canis lupus familiaris Major allergen Can f 1 protein
Target	Major allergen Can f 1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O18873
Expiration Date	12 months from date of receipt.