

Product Datasheet

Alpha-mammal toxin AaH2 antibody (orb862296)

Description	Rabbit polyclonal antibody to Alpha-mammal toxin AaH2. All of our unique antibodies are supplied with...
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Androctonus australis Alpha-mammal toxin AaH2 protein
Target	Alpha-mammal toxin AaH2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P01484
Expiration Date	12 months from date of receipt.