

Product Datasheet

cyp c 1.02 antibody (orb862375)

Description	Rabbit polyclonal antibody to cyp c 1.02.All of our unique antibodies are supplied with the followin...
Species/Host	Rabbit
Reactivity	Fish
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Cyprinus carpio cyp c 1.02 protein
Target	cyp c 1.02
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8UUS2
Expiration Date	12 months from date of receipt.