

Product Datasheet

MB21D1 antibody (orb862473)

Description	Rabbit polyclonal antibody to MB21D1.All of our unique antibodies are supplied with the following: 2...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Homo sapiens CGAS protein
Target	CGAS
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8N884
Expiration Date	12 months from date of receipt.