

Product Datasheet

Collagen IV antibody (orb862767)

Description	Rabbit polyclonal antibody to Collagen IV. All of our unique antibodies are supplied with the followi...
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Mus musculus Col4a2 protein
Target	Col4a2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P08122
Expiration Date	12 months from date of receipt.