

Product Datasheet

Outer capsid glycoprotein VP7 antibody (orb862899)

Description	Rabbit polyclonal antibody to Outer capsid glycoprotein VP7.All of our unique antibodies are supplie...
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Rotavirus A Outer capsid glycoprotein VP7 protein
Target	Outer capsid glycoprotein VP7
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P04511
Expiration Date	12 months from date of receipt.