

Product Datasheet

SH3PXD2A antibody (orb863042)

Description	Rabbit polyclonal antibody to SH3PXD2A.All of our unique antibodies are supplied with the following:...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Homo sapiens SH3PXD2A protein
Target	SH3PXD2A
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q5TCZ1
Expiration Date	12 months from date of receipt.