

Product Datasheet

map-1 antibody (orb863224)

Description	Rabbit polyclonal antibody to map-1. All of our unique antibodies are supplied with the following: 20...
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Mycobacterium tuberculosis map-1 protein
Target	map-1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P9WK20
Expiration Date	12 months from date of receipt.