

Product Datasheet

PLA2G12A antibody (orb863496)

Description	Rabbit polyclonal antibody to PLA2G12A. All of our unique antibodies are supplied with the following:...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Homo sapiens PLA2G12A protein
Target	PLA2G12A
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9BZM1
Expiration Date	12 months from date of receipt.