

Product Datasheet

Toxin Tb2-II antibody (orb863649)

Description	Rabbit polyclonal antibody to Toxin Tb2-II. All of our unique antibodies are supplied with the follow...
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Tityus bahiensis Toxin Tb2-II protein
Target	Toxin Tb2-II
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P60276
Expiration Date	12 months from date of receipt.