

Product Datasheet

Wt1 antibody (orb863774)

Description	Rabbit polyclonal antibody to Wt1. All of our unique antibodies are supplied with the following: 200 ...
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Mus musculus Wt1 protein
Target	Wt1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P22561
Expiration Date	12 months from date of receipt.