

Product Datasheet

L1R antibody (orb863945)

Description	L1R antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Vaccinia virus Protein L1(L1R),partial (2-183aa)
Target	L1R
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P20540
Expiration Date	12 months from date of receipt.