

Product Datasheet

Tmsb10 antibody (orb864089)

Description	Rabbit polyclonal antibody to Tmsb10. All of our unique antibodies are supplied with the following: 2...
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Mus musculus Tmsb10 protein
Target	Tmsb10
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q6ZWY8
Expiration Date	12 months from date of receipt.