

Product Datasheet

CCL17 antibody (orb864141)

Description	Rabbit polyclonal antibody to CCL17. All of our unique antibodies are supplied with the following: 20...
Species/Host	Rabbit
Reactivity	Canine
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Canis familiaris CCL17 protein
Target	CCL17
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q95N01
Expiration Date	12 months from date of receipt.