

Product Datasheet

GUCY1A2 antibody (orb864156)

Description	Rabbit polyclonal antibody to GUCY1A2. All of our unique antibodies are supplied with the following: ...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Homo sapiens GUCY1A2 protein
Target	GUCY1A2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P33402
Expiration Date	12 months from date of receipt.