



Product Datasheet

Angiotensin-converting Enzyme 2 protein [Out of stock] (orb864236)

Description	Recombinant human Angiotensin-converting Enzyme 2 protein
Reactivity	Human
Tested Applications	Functional Studies, IA, WB
Preservatives	PBS 20% Glycerol
Form/Appearance	Liquid
Storage	-20° or -70°C
Tag	10xHis at the C-terminus
Note	For research use only
Protein Sequence	19-740 aa
Purity	> 95%
Source	HEK293 cells
MW	85 kDa
Uniprot ID	Q9BYF1
Expiration Date	6 months from date of receipt.



Angiotensin-converting Enzyme 2 protein



Angiotensin-converting Enzyme 2 protein