


Product Datasheet

SARS-CoV-2 trimeric soluble full-length Spike protein, Beta variant, S protein (orb864237)


Description

Recombinant human SARS-CoV-2 trimeric




SDS-PAGE (left)
Western blot (right)

SARS-CoV-2 trimeric soluble full-length ...



SARS-CoV-2 trimeric soluble full-length ...



SARS-CoV-2 trimeric soluble full-length ...

Reactivity

variant, S protein

Human

Tested Applications

Functional Studies, IA, WB

Immunogen

Unconjugated

Preservatives

PBS 20% Glycerol

Form/Appearance

Liquid

Storage

-20° or -70°C

Tag

8xHis at the C-terminus

Note

For research use only

Protein Sequence

16-1211

Purity

> 95%

Source

HEK293 cells

MW

137

Uniprot ID

[P0DTC2](#)

Expiration Date

6 months from date of receipt.