

Product Datasheet

SARS-CoV-2 trimeric soluble full-length Spike protein, B1.617 cluster, S protein (orb864239)

Description Recombinant human SARS-CoV-2 trimeric soluble full-length Spike protein, B1.617 cluster, S protein

Reactivity Human

Tested Applications Functional Studies, IA, WB

Immunogen Unconjugated

Preservatives PBS 20% Glycerol

Form/Appearance Liquid

Storage -20° or -70°C

Tag 8xHis at the C-terminus

Note For research use only

Protein Sequence 16-1211


Purity > 95%

Source HEK293 cells

MW 137

Uniprot ID **P0DTC2**

Expiration Date 6 months from date of receipt.



SARS-CoV-2 trimeric soluble full-length ...