

Product Datasheet

Human GLP-1 (glucagon-like peptide-1) ELISA kit (orb864483)

Description	Human GLP-1 (glucagon-like peptide-1) ELISA kit
Reactivity	Human
Tested Applications	ELISA
Storage	4°C
Note	For research use only
Application notes	Detection samples: 1) Body fluid: Serum, plasma, tissue homogenate and other biological fluids; 2)Non-body fluid: cell culture supernatant, cells, small molecule peptide NB: Please specify your choice of sample types to our team.
Expiration Date	Please enquire.