

## Product Datasheet

### Mouse Epiregulin ELISA Kit (orb864913)

|                            |                                                                                                                                                                                                                                       |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Mouse Epiregulin ELISA Kit                                                                                                                                                                                                            |
| <b>Reactivity</b>          | Mouse                                                                                                                                                                                                                                 |
| <b>Tested Applications</b> | ELISA                                                                                                                                                                                                                                 |
| <b>Storage</b>             | 4°C                                                                                                                                                                                                                                   |
| <b>Note</b>                | For research use only                                                                                                                                                                                                                 |
| <b>Application notes</b>   | Detection samples: 1) Body fluid: Serum, plasma, tissue homogenate and other biological fluids; 2)Non-body fluid: cell culture supernatant, cells, small molecule peptide NB: Please specify your choice of sample types to our team. |
| <b>Expiration Date</b>     | Please enquire.                                                                                                                                                                                                                       |