

Product Datasheet

Mouse TNFRSF13C ELISA Kit (orb865100)

Description	Mouse TNFRSF13C ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Storage	4°C
Note	For research use only
Application notes	Detection samples: 1) Body fluid: Serum, plasma, tissue homogenate and other biological fluids; 2)Non-body fluid: cell culture supernatant, cells, small molecule peptide NB: Please specify your choice of sample types to our team.
Expiration Date	Please enquire.