


Product Datasheet

Human FCRL2 ELISA Kit (orb865264)

Description	Human FCRL2 ELISA Kit
Reactivity	Human
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0, Immunogen sequence: L20-V401
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA).
Uniprot ID	Q96LA5
Sensitivity	50 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is 250-550 pg/ml.
Expiration Date	Please enquire.



Human
FCRL2
PicoKine
ELISA Kit
Standard ...