


Product Datasheet

Mouse FCRL5 ELISA Kit (orb865265)

Description	Mouse FCRL5 ELISA Kit
Reactivity	Mouse
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0, Immunogen sequence: Q27-A496
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA).
Uniprot ID	Q68SN8
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is 100-200 pg/ml.
Expiration Date	Please enquire.



Mouse
FCRL5
PicoKine
ELISA Kit
Standard ...