


Product Datasheet

Human BAFF/TNFSF13B ELISA Kit (orb865266)

Description	Human BAFF/TNFSF13B ELISA Kit
Reactivity	Human
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: A134-L285
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA)
Uniprot ID	Q9Y275
Sensitivity	2 pg/ml
Expiration Date	Please enquire.



Human
TNFSF13B
PicoKine
ELISA Kit
standa...