


Product Datasheet

Human APOA1/Apolipoprotein A-I ELISA Kit (orb865269)

Description	Human APOA1/Apolipoprotein A-I ELISA Kit
Reactivity	Human
Range	3.12 ng/ml - 200 ng/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: D25-Q267
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA) and urine
Uniprot ID	P02647
Sensitivity	50 pg/ml
Expiration Date	Please enquire.



Human
APOA1
PicoKine
ELISA Kit
standard ...