


Product Datasheet

Human TNFRSF13C/BAFFR ELISA Kit (orb865270)

Description	Human TNFRSF13C/BAFFR ELISA Kit
Reactivity	Human
Range	78 pg/ml - 5,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: S7-A71
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA)
Uniprot ID	Q96RJ3
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Human
TNFRSF13C
PicoKine
ELISA Kit
stand...