


Product Datasheet

Mouse CXCL4/PF4 ELISA Kit (orb865287)

Description	Mouse CXCL4/PF4 ELISA Kit
Reactivity	Mouse
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: V30-S105
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA)
Uniprot ID	Q9Z126
Sensitivity	1 pg/ml
Expiration Date	Please enquire.



Mouse PF4
PicoKine
ELISA Kit
standard
cu...