


Product Datasheet

Human CXCL13/BLC/BCA1 ELISA Kit (orb865292)

Description	Human CXCL13/BLC/BCA1 ELISA Kit
Reactivity	Human
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: V23-R94
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA) and saliva
Uniprot ID	O43927
Sensitivity	1 pg/ml
Expiration Date	Please enquire.



Human
CXCL13
PicoKine
ELISA Kit
standard...