


Product Datasheet

Human Chemerin/RARRES2 ELISA Kit (orb865301)

Description	Human Chemerin/RARRES2 ELISA Kit
Reactivity	Human
Range	78 pg/ml - 5,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: V17-S163
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA)
Uniprot ID	Q99969
Sensitivity	20 pg/ml
Expiration Date	Please enquire.



Standard curve of Human Chemerin ELISA K...