


Product Datasheet

Mouse Complement Factor H/CFH ELISA Kit (orb865305)

Description	Mouse Complement Factor H/CFH ELISA Kit
Reactivity	Mouse
Range	0.78 ng/ml - 50 ng/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: S875-V1252
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA)
Uniprot ID	P06909
Sensitivity	20 pg/ml
Expiration Date	Please enquire.



Mouse CFH
PicoKine
ELISA Kit
standard
cu...