


Product Datasheet

Human CCL24/Eotaxin-2 ELISA Kit (orb865308)

Description	Human CCL24/Eotaxin-2 ELISA Kit
Reactivity	Human
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: V27-C119
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA)
Uniprot ID	O00175
Sensitivity	2 pg/ml
Expiration Date	Please enquire.



Human
CCL24
PicoKine
ELISA Kit
standard ...