


Product Datasheet

Mouse CCL12/MCP5 ELISA Kit (orb865318)

Description	Mouse CCL12/MCP5 ELISA Kit
Reactivity	Mouse
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: G23-G104
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA)
Uniprot ID	Q62401
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Mouse
CCL12
PicoKine
ELISA Kit
standard ...