


Product Datasheet

Mouse CD40/TNFRSF5 ELISA Kit (orb865320)

Description	Mouse CD40/TNFRSF5 ELISA Kit
Reactivity	Mouse
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: V24-R193
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA)
Uniprot ID	P27512
Sensitivity	1 pg/ml
Expiration Date	Please enquire.



Mouse CD40
PicoKine
ELISA Kit
standard c...