


Product Datasheet

Human KIT/SCFR/CD117/C Kit ELISA Kit (orb865326)

Description	Human KIT/SCFR/CD117/C Kit ELISA Kit
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: Q26-T520
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA)
Uniprot ID	P10721
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Human KIT
PicoKine
ELISA Kit
standard
cu...