


Product Datasheet

Human ASAH2 ELISA Kit (orb865328)

Description	Human ASAH2 ELISA Kit
Reactivity	Human
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: CHO; Immunogen sequence: F99-I780
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA)
Uniprot ID	Q9NR71
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Human
ASAH2
PicoKine
ELISA Kit
standard ...