


Product Datasheet

Human CCL15/Mip 1 Delta ELISA Kit (orb865334)

Description	Human CCL15/Mip 1 Delta ELISA Kit
Reactivity	Human
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: Q22-I113
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA)
Uniprot ID	Q16663
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Human
CCL15
PicoKine
ELISA Kit
standard ...